

FIG. 1

PURIFICATION SCHEME

1. Purification Step

Plasma



Immuno Affinity Chromatography

IgG-eTTP Affi-Gel® HZ 200 ml



Elution of vWF protease with NaSCN 3M



Gelfiltration

Bio-Gel® P-6DG 200 ml



Protein G-Sepharose ® 10 ml



Eluate

FIG. 1

2. Purification Step

Pooled Eluate from Step 1



Immuno Affinity Chromatography

Therasorb® 50 ml



Dilution



Anion Exchange Chromatography

High Q Support® 5 ml



Affinity Chromatography

Lentil Lectin-Sepharose® 25 ml



Elution of vWF Protease with Methyl- α -D-Mannopyranoside 0.3M

FIG. 1

3. Purification Step

Pooled Eluate from Step 2



Affinity Chromatography

Heparin-Sepharose 5ml



Dialysis



Anion Exchange Chromatography

High Q Support® 5 ml



Gel Filtration

Sephacryl® S-300 HR 500 ml



Elution of vWF Protease with NaCl 0.5 M

FIG. 1

Optional: 4. Purification Step

Pooled Eluate from Step 3



Affinity Chromatography

Anti- α_2 -Macroglobulin-column 1ml



and/or optional



Affinity Chromatography

Anti-Clusterin column

gctgcaggcgccatcctacacctggagctgctggtggccgtgggccccgatgtc
 A . A G G I L H L E L L V A V G P D V
 ttccaggctcaccaggaggacacagagcgctatgtgctcacccaacctcaacatcggggca
 F Q A H Q E D T E R Y V L T N L N I G A
 gaactgcttcgggacccgtccctgggggctcagtttcgggtgcacctgggtgaagatgggtcat
 E L L R D P S L G A Q F R V H L V K M V I
 tctgacagagcctgaggggtgccccaaatatcacagccaacctcacctcgctccctgctgagc
 L T E P E G A P N I T A N L T S S L L S
 gtctgtgggtggagccagaccatcaaccctgaggacgacacggatcctggccatgctgac
 V C G W S Q T I N P E D D T D P G H A D
 ctggctcctctatatcactaggtttgacctggagttgcctgatggtaaccggcaggtgcgg
 L V L Y I T R F D L E L P D G N R Q V R
 ggcgtcacccagctgggcggtgcctgctccccaacctggagctgcctcattaccgaggac
 G V T Q L G G A C S P T W S C L I T E D
 actggcttcgacctgggagtcaccatt
 T G F D L G V T I

147 AA

Fig. 2

147
 148
 149
 150
 151
 152
 153
 154
 155
 156
 157
 158
 159
 160
 161
 162
 163
 164
 165
 166
 167
 168
 169
 170
 171
 172
 173
 174
 175
 176
 177
 178
 179
 180
 181
 182
 183
 184
 185
 186
 187
 188
 189
 190
 191
 192
 193
 194
 195
 196
 197
 198
 199
 200